



**MEDICAL**

**DermaBlade<sup>®</sup>**  
**ASK FOR IT BY NAME**

- Flexible one-piece shave biopsy instrument
- Designed by a dermatologist for cutaneous surgery
- Smooth, clean excisions with less tissue trauma
- Individually packaged and sterilized

**Use for:**

- Removal of skin protuberances
- Saucerization
- Levelization
- Shave biopsy at any level of dermis

**Target Market**

Dermatologists  
 General Practitioners  
 Specialists



**Safe - Sterile**  
**Ready to Use**

ITEM CODE	PRODUCT DESCRIPTION	PACKING	SHIP UNIT	QTY PER SHIP UNIT	SHIP WT. lbs.	DIMENSIONS LxWxH in.
72-0001	DermaBlade <sup>®</sup> sterile	50/CT 8 CT/CASE	CS	400	3.50	11.50x10.50x8.50



**For more information:** call customer service 800-457-2222

**To order:** call your Medical Distributor

**To see the full line of Personna<sup>®</sup> medical products:** <http://medical.personnablades.com>

**Email:** [medicalmail@energizer.com](mailto:medicalmail@energizer.com)

**See video demonstration:** <http://dermablade.personnablades.com>

Personna<sup>®</sup>, DermaBlade<sup>®</sup> and MicroCoat<sup>®</sup> are registered trademarks of Energizer/Personna Industrial.

## Features

Surgical Instrument designed for shave biopsy.

## Flexibility

Sure grip side teeth/Three digit grip

MicroCoated<sup>®</sup> surgical blade,  
i.e., Persona Plus<sup>®</sup>

Controlled excision (shave biopsy)

## Safety

Sterile, individually packaged, with clear back

Disposable

Case size and quantity

Warning: Blade is extremely sharp. Handle with care.  
Used sharps are contaminated. Disease may occur from  
contact or injury. Dispose of in puncture resistant container.  
Refer to OSHA or local standards setting organizational  
guidelines for handling and disposal of sharps.

## Benefits

- Safer for nurse and practitioner
- Avoid cuts from contaminated blades
- No breaking of razor blades
- Professional appearance reduces patient anxiety
- Clean blade surface – no foreign material (oil)
  
- Use to remove surface protuberances or deep lesions (bends to accommodate required angle)
- Reduces ragged skin edges
  
- Easy to pick up
- Prevents slipping in gloved hand, regardless of blade position
- Provides stability and control
  
- Sharp, smooth, resistance-free cutting
- Quicker healing and reduced scarring
- Low profile instrument
  
- No sutures
- No needle holders, pick-ups, etc.
- Effective, but less costly set-up
  
- Easy to pick up from surgical tray
- No hemostat need to pick up blade
- Work place safety is improved (no razor blades to break, wrap, and use)
  
- Patient ready, sterility guaranteed
- Convenient, safe, easy to open
- Distinctive blue color increases product visibility
  
- Eliminates need to install or remove blade
- Single use to prevent cross contamination
- Enhanced user safety
  
- Sized to complement office needs
- Dispensing carton, fits on shelf